

ം ക

%

ZYMO RESEARCH The Beauty of Science is to Make Things Simple

Introducing foreign DNA to yeast cells is essential for two-hybrid system screening and many other genetic manipulation techniques. The transformation efficiency of yeast cells is typically magnitudes lower than that of E. coli due to the added challe -nges associated with permeating the fungal cell wall.

Explore the top tips and tools for maximizing your yeast transformations: https://lnkd.in/gi6JSFH9

#Yeast #LifeScience #Biotech #Biotechnology #Biology #ZymoResearch



How to Maximize Yeast Transformation Efficiency





03 Heat Shocking

⁰⁴ Plating Media

05 DNA Input

CLICK ON THE LINK TO DISCOVER TOP TIPS



805, Ashok Bhawan, 93 Nehru Place, New Delhi-110019 Tel.: (011) 41608700, 9810306768 | Email: info@everonlife.com | Web: www.everonlife.com